








GENESI versions with air dryer








30/10 GENESI M



Code: D-GE150-030-M

Description: Oil free compressor, FAD at 5 bar 120 l/min, tank capacity 30 l., with membrane air dryer free of maintenance and 0,01 micron filtration.

Model	 Dental chairs	 Air dryer	 Air Capacity	 FAD @ 5bar	 Motor	 Volt. Freq.	 Receiver
			m3/h l/min cfm	m3/h l/min cfm	HP KW	V Hz	l GAL
30/10 GENESI M	2	yes	- 168 6	- 120 4,2	1,5 1,1	230 50	30 7,92

Model	 Noise Level	 Duty cycle	 Max pressure	 Size	 Packaging	 Net weight	 Gross weight
	dB(A)		bar PSI	mm	mm	Kg lb	Kg lb
30/10 GENESI M	69	50%	7 101	570x430x620	572x432x622	37 81	39 85

