









PRIME VERSIONS WITH AIR DRYER






30/15 PRIME M



Code: PR-OF200-030-M

Description: Oil free compressors ,FAD at 5 bar 152 l/min, tank capacity 30 l, with internal surface treatment against corrosion, equipped with the new membrane air dryer able to warranty till -20

Model	 Dental chairs	 Air dryer	 Air Capacity	 FAD @ 5bar	 Motor	 Volt.	 Freq.	 Receiver
			m ³ /h l/min cfm	m ³ /h l/min cfm	HP KW	V	Hz	l GAL
30/15 PRIME M	2-3	yes	- 255 9,01	- 152 5,4	2 1,5	230	50	30 7,92

Model	 Noise Level	 Duty cycle	 Max pressure	 Size	 Net weight	 Gross weight
	dB(A)		bar PSI	mm	Kg lb	Kg lb
30/15 PRIME M	69	70%	7 101	680x400x650	53 117	54 119

